20913_Masajeador Facial Lisey

FRONT PART

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BACK PART

20913 - Face Massager Lisey User manual (Eng)

It is recommended to read this user manual carefully and keep it for future reference.

Product overview:

Especifications:

- Voltage: 3V==Battery CR2032(included)
- Vibration frequency: 8500RPM

Operation:

Remove the bottom compartment cover and remove the battery protector. If the battery is replaced, replace it with one of the same type, respecting the polarity.

Close the compartment and press the up/down switch to turn on the device.

In the upper (+) position the device will produce vibration and a positive current.
In the lower (-) position the device will produce vibration and a negative current.

Turn off the switch to the center position to turn off the device.

Warnings:

- The appliance is not intended for use by children or persons with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless such persons are given supervision or instructions necessary for the operation of the appliance. Keep away from children.
- Do not use this appliance near bathtubs, showers or water tanks.
- The use of this device may cause discomfort or adverse effects to people suffering from skin diseases or inflammations, people with sensitive skin and people with allergic constitution. Stop using it if you feel any discomfort.
- This device should not be used in areas where the skin is damaged.
- This device should not be used by people who have implants incompatible with this use.
- This device should not be used by persons under 14 years of age.
- Failure to follow instructions may cause adverse effects.
- Prolonged or repeated application to the same area could damage the skin.
- Do not override or remove the security feature inherent to the application.
- This device should not be used for the treatment of certain parts of the body, such as eyes or mucous membranes.

- Turn off the output of the device before placing or removing the electrodes on or from the skin.
- Do not use the device with flammable-based cosmetic products unless recommended by the manufacturer.
- Do not use the device near your eyes.
- Do not use the device on people with an implanted electronic device.
- Do not use the device near the heart.
- Do not expose it to rain or wet environments. Keep the device dry.
- Avoid falls that could cause serious damage.
- Do not attempt to repair the unit yourself, contact a specialist.
- To avoid risks, do not dispose of it inappropriately.
- Wipe it with a dry cloth.
- Do not install near heat sources.
- Do not insert any metal objects into the unit, there is a possibility of a short circuit.

This symbol both in the product and its user manual indicates that at the end of the devices electrical life, it should be recycled separately from your household waste. There are adequate means to dispose of this material for proper recycling. For more information, contact your local authorities.

MKTO CATAL IMPORTACIONES S.L. - ESB04267225
Ctra. Huércal-Overa, s/n - 04640, Puñal, Almería, Spain
productsafety@makitoe.es